ARJO

ARJO

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Contents

A Parts referred to in these instructions (Page 2)

B Introduction (Page 3)

C Care of Your Batteries (Page 4 to 5)


The vertical and horizontal lines printed in the margins adjacent to the text/illustrations in these instructions are for ARJO use only and should be disregarded by the reader.

Symbols used adjacent to the text in these instructions:-

-  **Danger** Means:- electrical hazard warning, failure to understand and obey this warning may result in electrical shock.
-  **Warning** Means:- failure to understand and obey this warning may result in injury to you or to others.
-  **Caution** Means:- failure to follow these instructions may cause damage to all or parts of the system or equipment.
-  **Note** Means:- this is important information for the correct use of this system or equipment.

A

Parts referred to in these Instructions

Arjo Lead/Acid Batteries

Battery Label


Read and thoroughly understand the operating instructions issued for the product.

Battery voltage DC and capacity in Amp hours.



Requires separate collection: this product should not be disposed of through the conventional waste disposal channels. Principle chemistry of the product.

Represents a recyclable product and principle chemistry of the product.

Name and address of the manufacturer.

24VDC • 4Ah

Manufactured for
ARJO Ltd.,
St. Catherine Street,
Gloucester GL1 2SL
U.K.
Tel. +44 (0)1452 500200

P1351_KKX12010.0 1sc.2

Fig. 1

These care instructions are intended for operators when using Arjo Lead/Acid batteries only and the instructions are designed to compliment the operating instructions issued for your ARJO lifter. When using ARJO lifters please ensure that you are familiar with these care instructions and all other operating instructions relating to ARJO lifters, chargers and ancillary equipment.

ARJO Lead/Acid battery packs are designed specifically for your ARJO lifter and are robustly constructed from high impact resistance plastic. Each pack contains valve regulated lead acid type (VRLA) batteries and require the minimum amount of care to maintain their working efficiency.

Your battery pack is an integral and necessary part of your ARJO lifter. The effective use and life expectancy of your battery will be greatly enhanced if the following protocol is respected.

With a fully charged battery pack it is possible to complete from between twenty and one hundred complete lifting cycles, with a patient, this will depend on which ARJO lifter is being used, the weight of the patient and the age and condition of the battery.

Recharging a battery pack before it reaches a low state of charge will prolong its life.

Lead acid batteries are classified as consumable items which means that the number of lifts from a fully charged battery pack will gradually reduce with age.

Your ARJO battery pack works at maximum efficiency at 20°C (68°F) (normal room temperature) at higher temperatures the battery will appear less efficient and take longer to charge.

For battery charging instructions always refer to the specific operating instructions for your lifter and any separate ARJO battery charging instructions (if applicable).

Some information contained in these instructions may become outdated, due to improvements made to this product in the future. If you have any questions regarding these instructions or your lifter, please contact ARJO or their approved distributor.

ARJO's policy is one of continuous development, and therefore reserves the right to change specifications without notice.

Battery Care

Battery Charging Protocol

- ! When first received charge the battery for a minimum of twenty four hours.
- ! To maximise the longevity of the battery and charger it is best practice to always disconnect the charger from the main electrical supply before inserting or removing the battery from the charger.
- ! Always thoroughly read and understand the lifter operating instructions and any separate battery charging instructions before charging the battery pack.
- ! Once fully charged and in use the battery should be recharged at either:-
 - A. The end of every work shift or
 - B. When the low battery warning light illuminates, or
 - C. Daily - as a **minimum**.

Whichever occurs first.

- ! Always charge the battery pack for at least eight hours.
- ! Every seven days leave the battery on charge for a minimum of fifteen hours to maximise longevity of battery.
- **Note:** *The battery pack may be left connected to the charger unit when it is fully charged without being damaged by overcharging, there is no maximum time limit. This will also ensure the battery is kept fully charged.*

A simple battery management protocol suggestion is to allocate two or three battery packs for each lifter, identify each battery pack by marking with one of the different stickers provided with this booklet, illustrated below:-


Fig. 2

Alternatively if you wish to use your own rotational system, you could use the blank adhesive stickers and mark accordingly.

- ! Do not continue to operate your ARJO lifter with the low battery warning light illuminating, unless in an emergency, as this will greatly reduce the life of the battery pack.
- ! Always provide a fully charged battery pack for the following work shift regardless of the number of lifts that have been performed, this will ensure maximum battery life and optimal lift performance.

The simple icon self adhesive label provided with this booklet may be removed and positioned at a convenient location, close to the battery charger for example, as a reminder to staff about the importance of regular and sequential charging protocol.

Handling of Batteries

- ! Be careful not to drop batteries as they could cause personal injury or become damaged.
- ! Do not drop the battery onto a hard surface this can damage the cells and components inside, even if the plastic casing looks undamaged dropping the battery may result in poor or no performance at all.

Storage of batteries

- ! If storing a battery pack for any length of time it is very important to fully recharge the battery before storage, and charge again before use.
- ! To prolong shelf life without intermittent charging store battery packs cool, ideally at 10°C (50°F) or less, but not below -40°C (-40°F).
- ! Keep the battery packs away from direct sunlight and heat sources e.g. radiators, elevated ambient temperatures above 30°C (86°F) will shorten service life of batteries.

(Continued...)

Disposal of Batteries

Always dispose of batteries in a responsible manner, to protect the environment contact your local authority for advice and disposal provisions for lead acid batteries.

General Safety Practices For battery packs

Only use the charger unit supplied with the lifter.

Do **NOT** charge battery packs in a sealed container.

Do **NOT** place batteries near, or dispose of, in a fire.

Do **NOT** place battery packs in close proximity to objects that produce sparks or flames.

Do **NOT** short circuit a battery.

Do **NOT** store batteries at temperatures in excess of 40°C (140°F).

Do **NOT** charge or store battery packs in the presence of inflammable anaesthetic gases.

Do **NOT** crush, puncture, open, dismantle or otherwise mechanically interfere with batteries.

Should the battery casing become cracked, and electrolyte come into contact with skin or clothing, wash immediately with water.

If the electrolyte contacts the eyes, wash immediately with copious amounts of water, and seek medical attention.

When disposing of battery packs and/or batteries, contact the appropriate local authority for advice.

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Contents

A Parts referred to in these instructions (Page 42)

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
Symbols used adjacent to the text in these instructions:

-  **Danger** Electrical hazard warning. Failure to understand and obey this warning may result in electrical shock.
-  **Warning** Failure to understand and obey this warning may result in injury to you or to others.
-  **Caution** Failure to follow these instructions may cause damage to one or all parts of the system or equipment.
-  **Note** This information is important for the correct use of this system or equipment.

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


Fig. 1

These care instructions are intended for operators when using ARJO Lead/Acid batteries only. They are designed to compliment the operating instructions issued for the ARJO lift. When using ARJO lifts, please ensure that you are familiar with these care instructions and all other operating instructions relating to ARJO lifts, chargers and ancillary equipment.

ARJO Lead/Acid battery packs are designed specifically for the ARJO lift and are robustly constructed from high impact resistance plastic. Each pack contains valve regulated lead acid type (VRLA) batteries which require the minimum amount of care to maintain their working efficiency.

The battery pack is an integral and necessary part of the ARJO Lift. The effective use and life expectancy of the battery will be greatly enhanced if the following protocol is respected.

With a fully charged battery pack, it is possible to perform between twenty and one hundred complete lifting cycles with a patient. The number of cycles will depend on which ARJO lift is being used, the weight of the patient and the age and condition of the battery.

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Lead acid batteries are classified as consumable items, which means that the number of lifts from a fully charged battery pack will gradually reduce with age.

The ARJO battery pack works at maximum efficiency at 20°C (68°F) (normal room temperature). At higher temperatures, the battery will be less efficient and take longer to charge.

To charge the battery, always refer to the specific operating instructions for the lift and any separate ARJO battery charging instructions (if applicable).

Some information contained in these instructions may become outdated, due to improvements made to this product in the future. If you have any questions regarding these instructions or the lift, please contact ARJO or their approved distributor.

ARJO's policy is one of continuous development and therefore, reserves the right to change specifications without notice.

Battery Care


Battery Charging Protocol

- ! When first received, charge the battery for a minimum of twenty-four hours.
- ! To maximize the longevity of the battery and charger, it is best practise to always thoroughly read and understand the lift operating instructions and any separate battery charging instructions before charging the battery pack.
- ! Once fully charged and in use, the battery should be recharged at either:
 - A. the end of every work shift, or
 - B. when the low battery warning light illuminates, or
 - C. daily - **as a minimum,**

whichever occurs first.

- ! Always charge the battery pack for at least eight hours.
- ! ARJO recommends that the battery be charged every 7 days for 15 hours. This will maximize the longevity of the battery.
- **Note:** *The battery pack may be left connected to the charger unit when it is fully charged, without being damaged by overcharging. There is no maximum time limit. This will also ensure that the battery is kept fully charged.*

A simple battery management protocol suggestion is to allocate two or three battery packs for each lift, and identify each battery pack by marking it with one of the different stickers provided with this booklet, as illustrated below:

Alternatively, if you wish to use your own rotational system, you could use the blank adhesive stickers and mark accordingly.

- ! Do not continue to operate the ARJO lift with the low battery warning light is illuminated, unless in an emergency, as this will greatly reduce the life of the battery pack.
- ! Always provide a fully charged battery pack for the following work shift, regardless of the number of lifts that have been performed. This will ensure the maximum battery life and optimal lift performance.

The self-adhesive icon label provided with this booklet, may be removed and positioned at a convenient location, i.e. close to the battery charger, as a reminder to staff about the importance of regular and sequential charging protocol.

Handling of Batteries

- ! Be careful not to drop batteries as they could cause personal injury or become damaged.
- ! Do not drop the battery onto a hard surface, as this can damage the cells and components inside. Even if the plastic casing looks undamaged, dropping the battery may result in poor performance or no performance at all.

Storage of batteries

- ! When stored, a battery will slowly self-discharge. Keep the stored battery packs in a dry environment at 20°C (68°F).
- ! To prolong shelf life without intermittent charging, store battery packs in a cool environment, ideally at 10°C (50°F) or less but not below -40°C (-40°F).
- ! Keep the battery packs away from direct sunlight and heat sources, i.e. radiators. Ambient temperatures above 30°C will shorten the service life of batteries.

Disposal of Batteries

Always dispose of batteries in a responsible manner. To protect the environment, contact your local authority for advice and disposal provisions for lead acid batteries.

General Safety Practices for battery packs

Only use the charger unit supplied with the lift.

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Do **NOT** place battery packs close to objects that produce sparks or flames.

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Do **NOT** charge or store battery packs in the presence of inflammable anesthetic gases.

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If the battery casing becomes cracked, and electrolyte comes into contact with skin or clothing, wash immediately with water.

If the electrolyte contacts the eyes, wash them immediately with copious amounts of water, and seek medical attention.

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

ARJO strongly advise and warn that only ARJO Company Designated Parts, which are designed for the purpose, should be used on equipment and other appliances supplied by ARJO, to avoid injuries attributable to the use of inadequate parts.

The ARJO Company's Conditions of sale make specific provision confirming no liability in such circumstances. ARJO's policy is one of continuous development, and therefore reserve the right to change specifications without notice.

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We strongly advise and warn that only Company Designated Parts, which are designed for the purpose, should be used on equipment and other appliances supplied by the Company, to avoid injuries attributable to the use of inadequate parts. The Company's Conditions of sale make specific provision confirming no liability in such circumstances.

Our policy is one of continuous development, and we therefore reserve the right to change specifications without notice.

AUSTRALIA
ARJO Hospital Equipment Pty. Ltd.,
154 Lytton Road
Bulimba
Brisbane
QLD, 4171
Tel. (61) 7 3395 6311
Fax. (61) 7 3395 6712

BELGIQUE / BELGIË
ARJO Hospital Equipment nv/sa
Twaalf Apostelenstraat 14
9051 Sint-Denijs-Westrem (Gent)
Tel. 09 242 86 42
Fax 09 242 86 43
E-Mail : info@ arjo.be

CANADA
ARJO Canada Inc.
1575 South Gateway Road
Unit C
Mississauga, ONTARIO
L4W 5J1
Tel. 1 800 665 4831
Fax. 1 800 309 1116

ČESKÁ REPUBLIKA
ARJO Hospital Equipment s.r.o.
Strmá 35
CZ-616 00 BRNO
Tel. 05 754 252
Fax. 05 412 13550

DANMARK
Getinge-ARJO A/S
Firkovvej 23
DK-2800 KGS Lyngby
Tel. 45 93 27 27
Fax. 45 93 41 20

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ARJO-Sic AG
ARJO Sverige AB
ARJO, Inc.
ARJO Limited • GLOUCESTER • U.K.