Mini L.E.D.

User's manual • Manuel d'utilisation
Manual de instrucciones • Bedienungsanleitung
Istruzioni per l'uso • Handleiding

ENGLISH	3
FRANCAIS	11
ESPAÑOL	19
DEUTSCH	27
ITALIANO	35
NEDERLANDS	43

CONTENTS

1	SAFETY INSTRUCTIONS	4
2	PRESENTATION	5
3	INSTALLATION	5
4	PRODUCT DESCRIPTION	6
5	MENUS	6
6	EVERYDAY USE	7
7	CARE	7
8	TROUBLESHOOTING	8
9	TECHNICAL SERVICE	8
10	SPECIFICATIONS	8
11	SYMBOLS AND ABBREVIATIONS	9

1 - SAFETY AND INSTRUCTION

Warning: Before connecting the Mini L.E.D., please read the safety instructions carefully. The manufacturer's liability for the Mini L.E.D. is applicable only if it is used as directed and if the safety instructions are properly complied with. Safety instructions mean a good prior knowledge of dentistry and photo-polymerization as well as understanding of the specific guidelines given in this User's Manual

- 1) This unit must be used in accordance with the instructions contained in this manual. We do not accept liability in the case of non compliance with the instructions or if the unit is used for any purpose other than that recommended.
- 2) Before connecting the base station to the mains, check that the mains voltage is set between the minimum and the maximum as indicated on the label on the mains adapter. A voltage outside the range stated could damage and/or injure the operator and/or the patient.
- 3) The light rays produced by the unit can be harmful and should never be directed towards the eyes even if adequate protective glasses are worn. The light is only intended for treatment inside the mouth.
- 4) Condensation inside an electrical appliance can be dangerous. If the unit needs to be transported from a cold to a warm area, it should not be used immediately but should be left long enough to reach the new room temperature.
- 5) Never replace the battery while the unit is plugged into the mains or before it has completely cooled off after use.
- 6) Do not use the unit near flammable anesthetic liquid, oxygen or nitrous oxide.
- 7) Unless otherwise specified in this manual, Do not introduce or try to introduce any metal objects into the unit or the base station. This could cause an electric shock.
- 8) We do not accept any liability for damage if accessories or replacement parts other than those issued by the manufacturer are used. In particular, the use of any other light guide, mains adapter or battery could be dangerous for the patient/user and/or for the correct operation of the instrument.
- 9) In the case of a fault, unplug the unit immediately and ensure that no other person can use it before it has been checked by the manufacturer or the manufacturer's authorized dealer. A fault could be due to incorrect use or physical damage.
- 10) Do not use near to any heat source. Use of solvent or flammable products on this unit could cause damage.
- 11) The use of detergents on the unit for cleaning or any other purposes may cause a short circuit.
- 12) If repairs are necessary, the unit should only be opened by the manufacturer's authorized dealer.
- 13) The **Mini L.E.D.** must not be used on any person with a heart stimulator or pacemaker.

- 14) Any unit used for polymerization, including the Mini L.E.D., must not be used on, or by, any person who suffers, or has suffered in the past, from photo-biological reactions (including urticaria solaris or protoporphyrie erythropetique), or any person receiving medical treatment and using photosensitive medication (including methoxsalenes or dymethylchlorotetracycline).
- 15) Any person who has previously suffered from a retina or crystalline condition or who has undergone eye surgery, in particular cataract treatment, should visit an eye specialist before using the unit or before undergoing treatment with the unit. Even with the patient's agreement, great care should be taken because the intensity of the light is potentially dangerous. It is highly recommended to use protective ultra-violet and light-filtering glasses at all times.
- 16) Any variation in the voltage of the power supply or in the electromagnetic field outside the normal safety regulations could cause the unit to come on automatically, or interfere with its operation.
- 17) Electromagnetic interference: the unit has been tested and approved for all the current standards (EN 60 601-1-2) both with regard to emission (electrical field rays and supply voltage interference) and immunity (electrical fields rays, intermittent bursts on the supply, electrostatic discharges and supply shock waves).
- 18) If required, your dealer can supply all spare parts, in particular the battery, which should be returned to them.
- 19) For use only by a dental professional.

2 - PRESENTATION

The Mini L.E.D. produces a visible blue light between a wavelength of 420 and 480nm for the photo-polymerization of dental materials. After checking the packaging for damage (please retain the packaging), you will find the following:

- 1- The Mini L.E.D. handpiece incorporating the Light Emitting Diode.
- 2- The base station for recharging the unit incorporating a radiometer.
- 3- The mains adapter and mains socket.
- 4- The Lithium Ion battery.
- 5- A curved universal Light Guide of Ø7,5mm which can be sterilized.
- 6- A deflector
- 7- Booster tip (optional)
- 8- User's manual.

3 - INSTALLATION

- 1) Check that environmental conditions are satisfactory (i.e. temperature between 15°C and 30°C and humidity between 10% and 80%).
- 2) Prior to each use, sterilize the light guide and disinfect both base station and handpiece (see also the chapter on care).
- 3) Screw the battery onto the handpiece, remove the protection cap and place the sterilized light guide in the handpiece. Ensure that the light guide is inserted correctly. There is a friction contact with a click.
- 4) Position the base station on a sturdy surface. Plug in the mains adapter and ensure that the plug is suitable for the voltage of your power supply. The green light will flash twice.
- 5) Place the handpiece in the base station for a minimum of 2 hours 30 minutes in order to fully charge the battery.
- 6) The green light starts to flash indicating that the battery is being charged. As soon as the battery is fully charged, the light stops flashing and stays lit.

4 - PRODUCT DESCRIPTION

Now that the **Mini L.E.D.** is fully charged, you need to familiarize yourself with its buttons and menus:

Handpiece (1):

LED indicators:

The handpiece has four LEDs:

- one status LED close to the light guide indicating the following:

Color	Mode
Green	Normal mode
Red	Battery low
Red flashing	Thermal protection

- three LEDs in line showing the selected curing mode:

Color	Mode
Red	Fast Curing
Orange	Pulse curing
Yellow	Step curing

Buttons:

The handpiece has two buttons:

- one navigation button which is used to switch from one menu to another
- One on/off light button (trigger) close to the light guide.

Base station (2):

Radiometer:

A built-in radiometer, associated with a LED indicator allows the operator to check accurately the curing unit's performance. The aperture is located on the front part of the base station.

To proceed:

- Set the unit in Fast Curing mode.
- Insert the standard light guide into the hand piece after checking no VLC material remains on the light guide.
- Place the end of the light guide in the aperture and activate the **Mini L.E.D.**

LED indicator:

Color	Mode
Green	Power over 800mW/cm ²
Red	See troubleshooting chapter

Charge:

When the handpiece is inserted in the base station, another LED displays the charge status for the battery:

Color	Mode
Green flashing	Charging
Green	Battery fully charged

5 - MENUS

The **Mini L.E.D.** has three different menus. To select the menu you want to use, press the navigation button.

Fast Curing Menu: when the handpiece is removed from the base station, the unit will be set on fast curing program for 10 seconds. The red light indicates your selection.

In this mode the power provided by your unit is approx.:

- 1100mW/cm² (± 10%) with the Ø7.5mm standard light guide.
- 2000mW/cm² (± 10%) with the Ø5.5mm optional booster Tip.

Pulse Curing Menu: selected by the navigation button. This menu gives full power in a pulsation mode with an emission of ten successive flashes including a rest period of 250ms between the flashes. The orange light indicates your selection.

In this mode the unit provides the same power described in the fast curing menu.

Step Curing Menu: Selected by the navigation button. If you wish to carry out "step curing" (curing similar to a step menu in halogen unit), you will have:

- soft start for 10 seconds from 0 to 1100mW/cm² (or 0 to 2000mW/cm² with the optional booster Tip)
- Full power during 10 seconds.

6 - EVERYDAY USE

- 1) The **Mini L.E.D.** is normally positioned on the base station. The first time you pick it up it will automatically be activated in the fast curing menu. For further details, please refer to chapter 5.
- 2) When the preferred program has been selected, your **Mini L.E.D.** is ready for use.
- 3) Place the light guide as close as possible to the surface of the VLC material without touching it, as this could cause damage to the light guide and decrease its efficiency.
- 4) Press the trigger button briefly in order to set the curing cycle. The setting will be confirmed by an audible signal (beep).
- 5) At the end of the cycle, you will hear another "beep". However, you can interrupt the cycle at any time simply by pressing the trigger button.
- 6) After three minutes, if the unit is not used, it goes into stand-by mode and all the lights go out. Press any button to reactivate the unit.

7 - CARE

Warning: For everyday care of the **Mini L.E.D.**, ensure that the handpiece is not on the base station before you remove the light guide. Unplug the **Mini L.E.D.** before using a disinfectant. Clean and sterilize the light guide before each patient.

- 1) The **Mini L.E.D.** and accessories can be disinfected using normal surface disinfectant, for example, a hydro-alcohol solution. Ensure that the **Mini L.E.D.** is unplugged before disinfecting and do not use abrasive products on the casing.
- 2) After each use, check that no trace of VLC material remains on the light guide. Any remaining VLC material must be removed immediately and the surface of the light guide checked, ensuring that it has not been altered in any way. If it is damaged, replace the light guide as the power could be reduced by up to 30%.
- 3) Under normal conditions, the power of the LED source does not vary if the battery is correctly charged and consequently the power level does not need to be checked as with ordinary curing units. However, if in doubt, use the built-in radiometer described in chapter 4.
- 4) The battery is located at the end of the handpiece. It can be removed by unscrewing.
- 5) Regularly check that the LED is clean. It is to be found inside the end of the handpiece and the light guide is slotted over it. If it is not clean, use a special cleaning solution developed for this purpose.

8 - TROUBLESHOOTING

If a problem appears, before calling your dealer's after sales service:

- 1) Check that the base station is correctly connected to the mains to ensure that the battery recharges normally. If the LED is not lit up in spite of correct connection, contact the after sales service.
- 2) Check that the battery indicator shows that it is fully charged before pressing the trigger. If the LED is red, an automatic safety mechanism will stop the unit from being activated. In this case, recharge your battery by replacing the handpiece on the base station, or just replace with a second fully charged battery.
- 3) Intensive use of the unit may possibly cause a rise in temperature inside the unit. In this occurs, an automatic safety mechanism will stop the **Mini L.E.D.**, and the status indicator will flash red. Leave the unit for a few minutes to cool down and then carry on as normal.

9 - TECHNICAL SERVICE

Your dealer's after sales service is at your disposal for help with any technical problem.

10 - SPECIFICATIONS


Unit

Model:	Mini L.E.D.
Weight:	190g
Size:	Ø23 x 200mm
Classification:	Ordinary Type B Continuous service IPX0

Power supply adapter

Input voltage:	110 - 240V AC
Frequency:	47 - 63 Hz
Output voltage:	12V
Output current:	1.25A
Classification:	II IP 40

Base station

Input voltage:	12VDC
Protection:	Fuse 2A
Classification:	Continuous service IPX 0

Battery

Type:	Lithium-Ion
Size:	90 x Ø21mm
Capacity:	2000mAh



Optical specifications

Wavelength:	420-480nm
Intensity (Ø7.5mm):	1100 mW/cm ²
Power:	450-500 mW

Homologation

CE
IEC 60 601-1
60 601-1-2

11 - SYMBOLS AND ABBREVIATIONS

SYMBOL	DEFINITION
	Alternating current
	Continuous current
	Warning, please refer to the accompanying documents
	Type B
	Class II

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