


■ Becton Dickinson Clay Adams* Dynac* and Dynac III Centrifuges

Compact, versatile centrifuge for clinical or industrial applications



Fits easily inside refrigerator or coldroom for low-temperature use.

- Variable speed control from 500 to 3600rpm/2224xG
- For volumes up to 100mL
- Built-in tachometer
- Timed-run (up to 30 min.) and continuous-run modes
- Electronic brake
- Positive rotor locking secures rotor on shaft to prevent free spinning
- Safety latch secures lid until rotor comes to complete stop
- Heavy-gauge aluminum guard bowl, baked enamel housing
- Transparent, shatterproof Lexan* cover

Ordering Information: For 110/120V 50/60Hz; accessory rotors available online.

Includes: Centrifuge only.

Required Accessories: Rotor or rotor assembly (with rings and buckets), all sold separately.

Model	Dimensions	Mfr. No.	Cat. No.	Each
Dynac III	24L x 20 7/8 W x 19" H (61 x 52 x 48cm)	420104	B420104	3955.79

■ Fixed-Angle Rotors for Dynac* and Dynac III Centrifuges

Six-, 12-, and 24-place fixed-angle rotors accept tubes from 3 to 15mL when equipped with buckets (sold separately). For optimum performance, the use of heavy-duty glassware is recommended in these rotors. For use with Dynac, Dynac II, and Dynac III Centrifuges.



Rotor System	BD No.	Cat. No.	Each
6-place	420111	05-100-31	819.38
12-place	420113	05-100-33	827.11
24-place	420114	05-100-34	827.11


■ 50mL Capacity 4-Place Horizontal Rotor for Dynac* and Dynac III Centrifuges

4 x 50mL capacity

- For use with: Dynac, Dynac II, and Dynac III Centrifuges
- Maximum speed/force: 2790rpm/1425xG in Dynac; 2690rpm/1397xG in Dynac III centrifuge

Includes: Rotor body and trunnion rings.

Required Accessories: Buckets (05-100-36, sold separately).



Capacity	BD No.	Cat. No.	Each
4 x 50mL	420108	05-100-28	827.11

■ 8-Place Horizontal Rotor for Dynac* and Dynac III Centrifuges

8 x 15mL capacity

- Rotor is supplied with trunnion rings to accommodate tubes from 3 to 15mL (including Vacutainer* tubes) when equipped with 05-100-37, -38, and -39 buckets (sold separately)
- For use with Dynac, Dynac II, and Dynac III Centrifuges.



BD No.	Cat. No.	Each
420109	05-100-29	827.11

■ Buckets for Fixed-Angle Rotors for Dynac* and Dynac III Centrifuges

For Tubes	O.D. x L	BD No.	Cat. No.	Each
12 to 15mL	17.6 x 118mm	420901	05-100-37	31.30
5 to 15mL	17.6 x 100mm	420902	05-100-38	43.87
3mL	13.8 x 75mm	420904	05-100-49	58.05