

Nor-Lake* Scientific Premier Glass-Door Freezers new

Choice of two capacities for laboratory and pharmacy storage



- Adjustable microprocessor control from -15°C to -25°C (factory preset at -20°C)
- Large digital display, accuracy to within 1°C
- High/low visual alarm
- Automatic defrost
- Top-mounted refrigeration system
- Thermostatic expansion valve meets varying load conditions
- Foamed-in-place, high-density polyurethane foam insulation
- Full glass doors with key locks, magnetic gasketing
- Air-cooled condensing unit
- Automatic condensate evaporator
- Four casters, two locking
- Door-opening perimeter heater wires
- Epoxy-coated shelves
- Interior fluorescent light
- White exterior and interior; galvanized steel top and bottom
- For 115/208-230V 60Hz

Compliance: UL and C-UL listed.

Note: Options such as audible high/low alarm with remote contacts, legs instead of casters, stainless-steel exterior/interior and more are available. Contact your Fisher representative for more information.

Warranty: 18-month parts and labor; 5-year compressor (US and Canada); 18-month parts (international).

Capacity	22.5 cu. ft. (637L)	48.8 cu. ft. (1382L)
No. of Doors	1	2
No. of Shelves	3	6
Casters	4	4
Exterior H x W x D	81 ¹ / ₂ x 27 ¹ / ₂ x 34 ⁷ / ₈ " (207 x 70 x 88cm)	83 ¹ / ₄ x 55 ¹ / ₂ x 34 ⁷ / ₈ " (212 x 140 x 88cm)
Mfr. No.	NSXF211WWG/8	NSXF482WWG/8
Cat. No.	22-650-631	22-650-632
Each	6722.00	9931.00

Nor-Lake* Scientific Premier Low-Temperature Freezers new

Adjustable temperatures from -30° to -20°C for safe, reliable storage



- 24 or 52 cu. ft. capacity
- Solid door(s)
- Digital LED microprocessor control with high/low visual alarm
- Foamed-in-place CFC-free high-density polyurethane foam
- CFC-free refrigerant
- Thermostatic expansion valve to meet varying load conditions
- Automatic defrost
- Centered key door locks
- Full-length door pull handles
- Magnetic door gaskets
- Epoxy-coated wire shelving
- Interior light
- Air-cooled condensing unit
- Automatic condensate evaporation
- White exterior, aluminum interior
- For 115/208-230V 60Hz

Compliance: UL listed.

Note: Optional features such as high/low audible alarm with remote contacts, temperature chart recorder, digital monitor thermometer with alarm, stainless-steel drawers, legs instead of casters, and more are available. Contact your Fisher representative for details.

Warranty: 18-month parts and labor; 5 years compressor (US and Canada); 18-month parts (International).

Capacity	24 cu. ft. (680L)	52 cu. ft. (1472L)
No. of Door(s)	1 solid	2 solid
No. of Shelves	3	6
Exterior H x W x D	81 ¹ / ₂ x 27 ¹ / ₂ x 34 ⁷ / ₈ " (207 x 70 x 88cm)	83 ¹ / ₄ x 55 x 34 ⁷ / ₈ " (211 x 140 x 88cm)
Mfr. No.	NSXF211WWW/8	NSXF482WWW/8
Cat. No.	22-650-629	22-650-630
Each	6205.00	8896.00