

■ Nor-Lake* Scientific Undercounter Refrigerators new

Versatile, energy-efficient refrigerators with choice of capacities and defrost modes

These space-saving refrigerators can be placed under the counter or on the benchtop; 22-650-600 can be stacked two high (with optional stacking kit) or left freestanding. Can be built in on three sides.

- Adjustable temperature control from 2° to 8°C (factory preset at 4°C)
- CFC-free refrigerant
- Foamed-in-place CFC-free high-density polyurethane foam insulation
- Pull door handle
- Key door locks
- Magnetic door gasketing
- Gray epoxy-coated shelving
- Base-mounted refrigeration system
- Air-cooled condensing unit
- Four 4" or two leveling legs, depending on model
- White enamel front, door, top, and sides; aluminum interior
- Galvanized steel back and bottom

Compliance: UL/CUL listed; 5.5 cu. ft. model also CEC (California Energy Commission) rated and Energy Star* compliant.

Warranty: 18 months parts and labor, 5 years compressor (US and Canada); 18 months parts (International).



Note: Optional features such as digital temperature display, high/low alarm with remote contacts, stacking kit, extra shelves, casters instead of legs, temperature chart recorder and more are available. Contact your Fisher representative for details.

Capacity	4.2 cu. ft. (119L)	4.6 cu. ft. (130L)	5.1 cu. ft. (144L)	5.4 cu. ft. (153L)	5.5 cu. ft. (156L)
Defrost	Automatic (forced air)	Manual (cold wall)	Automatic (forced air)	Manual (cold wall)	Automatic (forced air)
No. of Shelves	1	1	1	1	3 wire, 1 glass, 3 on door
No. of Drawers	0	0	0	0	1
Legs	Four 4" legs	Four 4" legs	Four 4" legs	Four 4" legs	Two front leveling legs
Exterior H x W x D	34-35 1/4 x 27 1/2 x 23 7/8" (86-90 x 90 x 61cm)	34-35 1/4 x 27 1/2 x 23 7/8" (86-90 x 90 x 61cm)	34-35 1/4 x 27 1/2 x 25" (86-90 x 69 x 64cm)	34-35 1/4 x 27 1/2 x 25" (86-90 x 69 x 64cm)	33 1/2 x 23 x 23 1/2" (85 x 58 x 60cm)
Mfr. No.	NSLR041WMW/0	NSLR041WMW/0M	NSLR051WMW/0	NSLR051WMW/0M	LR061WWW/0
Cat. No.	22-650-640	22-650-641	22-650-638	22-650-639	22-650-600
Each	2745.00	2644.00	2796.00	2758.00	1100.00