

Autoclaves

Thermo Scientific Barnstead SterileMax* Tabletop Sterilizers

17% more chamber space than traditional
benchtop units



- Holds upto three 2L flasks simultaneously, or several smaller items, using included removable racks and trays
- Four preset cycles allow one-touch operation for liquids, wrapped or unwrapped goods, and packs
- Built-in water pump maintains consistent water supply and accurate steam
- Chamber is insulated for quick cycle start-up and low heat loss between cycles
- LCD offers continuous viewing of the sterilization and diagnostic programs—plus, it displays faults to indicate that sterilization parameters have not been reached in chamber
- Integral access port allows additional verification of temperature and pressures when needed
- Available with a thermal printer that provides a paper tape with validation of sterilizer performance, load-control identification, and quality assurance checks for each cycle
- Optional B/T Sure* Biological Indicators for monitoring sterilizer performance
- Included three chamber locks will not open unless chamber has depressurized, as well as a low-water indicator display

Certifications: UL and CUL listed; 220V models carry CE Mark.

Compliance: ASME approved and stamped.

Electrical Requirements	Mfr. No.	Cat. No.	Each
With Printer			
120V 50/60Hz, 12.5A, 1500w	ST75935	14-488-24	7672.00
220V 50/60Hz, 6.5A, 1500w	ST75930 33	14-490-2	7672.00
Without Printer			
120V 50/60Hz, 12.5A, 1500w	ST75925	14-490-1	6905.00
220V 50/60Hz, 6.5A, 1500w	ST75920 33	14-490-3	6905.00